

GalaFLEX® scaffold is a bioabsorbable, monofilament surgical scaffold, constructed of poly-4-hydroxybutyrate (P4HB) – an advanced, biologically derived polymer that was developed by scientists at MIT.








# Discover the Next Generation of Soft Tissue Regeneration in Breast Surgery

GalaFLEX® scaffold is a biologically derived surgical scaffold that provides immediate mechanical support to the repair site. Over the course of approximately 24 months, GalaFLEX® bioabsorbs and encourages rapid tissue integration into the macropores of the monofilament design, strengthening tissue and resulting in a neotissue plane that is 3-4 times stronger than native tissue.<sup>9,16,19</sup> GalaFLEX® scaffold is designed to support, repair, elevate and reinforce soft tissue in the breast during surgical procedures<sup>1,4</sup> such as:

- Reduction mammoplasty
- Mastopexy
- Breast revision surgery

**GalaFLEX® scaffold offers a unique combination of properties that are optimal for soft tissue support in both medically necessary and cosmetic breast procedures:**

-  **Biologically Derived:** Produced by a **safe biological fermentation** process, standard in pharmaceutical production.<sup>12,17</sup>
-  **Monofilament:** Designed to minimize risk of **infection and encourage** a natural healing response.<sup>12,17</sup>
-  **Strong:** Provides a lattice for new tissue ingrowth and regeneration resulting in tissue **3-4x stronger than native tissue.**<sup>2,17</sup>
-  **Bioabsorbable:** Naturally broken down to CO<sub>2</sub> and H<sub>2</sub>O, with **bioabsorption essentially complete by 18-24 months.**<sup>1,12,17</sup>



## Comparative Scaffold Characteristics

	GalaFLEX® 2,20	VICRYL™ mesh 9,10	TIGR™ 2,7,14	STRATTICE™ 19,21
Material	P4HB	PLGA	PGLATMC/ PLATMC	Porcine
Structure	Monofilament	Multifilament	Multifilament	Acellular Dermal Matrix
Absorption Time (Months)	18-24	3	24-36	Remodels
Primary Absorption Mechanism	Hydrolytic	Hydrolytic	Hydrolytic	Enzymatic Remodeling
Initial Scaffold Bursts Strength (kgf) <sup>2</sup>	22.5	28.6	19.0	65
Retained Scaffold Strength at 12 weeks	>70%	0%	50%	21%

**Disclaimer** The above discussion points are in the context of the general literature, and not indicative of results from a head-to-head study.

**Intended Use** GalaFLEX® scaffold is intended for use, as an adjunct to sutures, for the reinforcement and repair of soft tissue where weakness exists and where the addition of a reinforcing material is needed to obtain the desired surgical result in patients undergoing breast surgery. The GalaFLEX® scaffold is designed to be used in patients undergoing soft tissue repair and reinforcement in medically necessary breast surgery procedures where the existing soft tissue is deficient to support the surgical repair. Examples of such breast surgery applications include reduction mammoplasty and breast revision surgery to correct a medical condition. GalaFLEX® scaffold may also be used in cosmetic breast procedures.

**Consult the GalaFLEX® Instructions for Use for complete prescribing information; including its indications for use, warnings and precautions.**

### References

- 1 Data on file at Tephra.
- 2 Preclinical data on file at Tephra; results may not correlate to clinical performance in humans.
- 3 Engelsman, A. F., van der Mei, H. C., Ploeg, R. J., & Busscher, H. J. "The Phenomenon of Infection with Abdominal Wall Reconstruction," Biomaterials, vol. 28 no. 14, 2018, pp. 2314-2327.
- 4 Adams Jr, W. P., Baxter, R., Glicksman, C., Mast, B. A., Tantillo, M., & Van Natta, B. W. (2018). "The Use of Poly-4-Hydroxybutyrate (P4HB) Scaffold in the Ptotic Breast: a Multicenter Clinical Study," Aesthetic Surgery Journal, 38(5), 502-518.
- 5 Alloderm Regenerative Tissue Matrix, Lifecell Electronic IFU, <https://allergan-web-cdn-prod.azureedge.net/actavis/>.
- 6 Amid, P. K., "Classification of Biomaterials and Their Related Complications in Abdominal Wall Hernia Surgery," Hernia, vol. 1, no. 1, 1997, pp. 15-21.
- 7 "Surgical Mesh Electronic IFU," TIGR Surgical, <http://novusscientific.com>.
- 8 Corey R Deeken, PhD, et al. "Histologic and Biomechanical Evaluation of Crosslinked and Non-Crosslinked Biologic Meshes in a Porcine Model of Ventral Incisional Hernia Repair," 23 Mar, 2011, <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3782991/>
- 9 Deeken, Corey R., and Brent D. Matthews, "Characterization of the Mechanical Strength, Resorption Properties, and Histologic Characteristics of a Fully Absorbable Material (Poly-4-Hydroxybutyrate—PHASIX Mesh) in a Porcine Model of Hernia Repair," ISRN surgery, 2013.
- 10 "Wound Closure Manual," Ethicon, Inc. PDF file, 2005.
- 11 Gross, John E., et al. "An Evaluation of SERI Surgical Scaffold for Soft-Tissue Support and Repair in an Ovine Model of Two-Stage Breast Reconstruction," Plastic and Reconstructive Surgery, vol. 134, no. 5, 2014, pp. 700e-704e.
- 12 Guo, Kai & Martin, David, "Poly-4-hydroxybutyrate (P4HB) in Biomedical Applications and Tissue Engineering," 2015.
- 13 Halaweish, Inab, et al. "Novel In Vitro Model for Assessing Susceptibility of Synthetic Hernia Repair Meshes to Staphylococcus Aureus Infection Using Green Fluorescent Protein-Labeled Bacteria and Modern Imaging Techniques," Surgical Infections, vol. 11, no. 5, 2010, pp. 449-454.
- 14 Hjort, H., et al. "Three-Year Results From a Preclinical Implantation Study of a Long-Term Resorbable Surgical Mesh with Time-Dependent Mechanical Characteristics," Hernia, vol. 16, no. 2, 2012, pp. 191-197.
- 15 Klinge U, Junge K, Spellerberg B, Piroth C, Klosterhalfen B, Schumpelick V, "Do Multifilament Alloplastic Meshes Increase the Infection Rate? Analysis of the Polymeric Surface, the Bacteria Adherence, and the In Vivo Consequences in a Rat Model," J Biomed Mater Res, vol. 63, no. 6, 2002, pp. 765-771.
- 16 Martin, D P, et al. "Characterization of Poly-4-Hydroxybutyrate Mesh for Hernia Repair Applications," Journal of Surgical Research, vol. 184, no. 2, 2013, pp. 766-773.
- 17 Martin, David P., and Simon F. Williams, "Medical Applications of Poly-4-Hydroxybutyrate: a Strong Flexible Absorbable Biomaterial," Biochemical Engineering Journal, vol. 16, no. 2, 2003, pp. 97-105.
- 18 Roth, J. S., et al. "Prospective Evaluation of Poly-4-Hydroxybutyrate Mesh in CDC Class L/High-Risk Ventral and Incisional Hernia Repair: 18-Month Follow-Up," Surgical Endoscopy, vol. 32 no. 4, 2013, pp. 1929-1936.
- 19 Scott, J. R., Deeken, C. R., Martindale, R. G., Rosen, M. J. "Evaluation of a Fully Absorbable Poly-4-Hydroxybutyrate/Absorbable Barrier Composite Mesh in a Porcine Model of Ventral Hernia Repair," Surgical Endoscopy, vol. 30, no. 9, 2016, pp. 3691-3701.
- 20 Williams, Simon F., Martin, David P., Moses, Arikha C, "The History of GalaFLEX P4HB Scaffold," Aesthetic Surgery Journal, 2016, pp. S33-S42.
- 21 "Mesh Electronic IFU," Stratice Surgical, <https://allergan-web-cdn-prod.azureedge.net>.

### GalaFLEX®: Available Sizes and Shapes

Shape	Product Code	Size (cm)
	CE0103	2.5 x 7.6
	CE0206	5 x 15
	CE0208	5.0 x 20.0
	CE0408	10 x 20
	CE0608	15.0 x 20.0







Tephra, Inc.  
99 Hayden Avenue, Suite 360  
Lexington, MA 02421, USA.  
[www.galateasurgicalinternational.com](http://www.galateasurgicalinternational.com)  
[contact-international@galateasurgical.com](mailto:contact-international@galateasurgical.com)



Galatea Surgical © 2020. All rights reserved. 500032 Rev H



# Strengthens Tissue in Breast Surgery

-  Biologically Derived
-  Monofilament
-  Strong
-  Bioabsorbable





# Strength and Beauty


## Biologically Derived

- Proprietary fermentation process designed and optimized to provide a biocompatible product that when combined with all other features encourages the patient's natural healing response.<sup>2,12,17</sup>
- P4HB devices have been tested in pre-clinical and clinical studies to ensure safety and effectiveness.<sup>2,18,19</sup>
- More than 3 million patients worldwide have had P4HB devices implanted.<sup>1</sup>

## Monofilament

- Designed with an open pore knit pattern to encourage rapid tissue ingrowth and to reduce risk of infection.<sup>3,6</sup>
- It has been reported that monofilament materials have on average 60% less surface area than that of multifilament materials, which may improve the healing response.<sup>3,15</sup>
- With less surface area, monofilament scaffolds have fewer recesses that bacteria can use as a haven from the body's natural defense systems or antibiotic treatments.<sup>3,13</sup>


When comparing SEM images of Galatea Scaffolds and other resorbable materials, the open pores, smooth surface and monofilament structure of Galatea Scaffolds are clearly visible.



**GalaFLEX® Scaffold**  
Monofilament  
derived from P4HB  
SEM Photo, 20x



**TIGR™ Mesh**  
Multifilament  
SEM Photo, 20x



**VICRYL™ Mesh**  
Multifilament  
SEM Photo, 20x


# Inside and Out

## Strong

- Designed specifically for strength retention throughout the critical wound healing period.<sup>1,16</sup>
- Rapid tissue regeneration resulting in a new tissue plane approximately 3-4 times the strength of the native tissue as demonstrated in pre-clinical studies.<sup>9,19</sup>
- Maintains >70% of its strength at 12 weeks in vivo.<sup>2</sup>

## Long-Term Repair Strength in a Preclinical Model<sup>9</sup>

(per Deeken, Matthews et al.)



## Bioabsorbable

- Naturally bioabsorb, leaving behind only strong, healthy tissue to support the surgical outcome.<sup>2,12</sup>
- Gradually and predictably bioabsorbs over the course of approximately 18-24 months.<sup>12</sup>
- Eliminated from the body as carbon dioxide and water primarily by the process of hydrolysis.<sup>9,12</sup>
- No polymer metabolites remain after the degradation process is complete.<sup>2</sup>

## Before Implantation



GalaFLEX® scaffold is a macroporous, monofilament, bioabsorbable scaffold.<sup>2</sup>

## After Implantation

(Human Breast Tissue Specimens)



Tissue rapidly grows into the pores of the GalaFLEX® scaffold, and forms a well-vascularized tissue plane.<sup>2</sup>




The newly formed tissue is pliable and provides strength and support to the elevated tissue.<sup>16</sup>

# GalaFLEX® encourages new tissue ingrowth and regeneration

- Provides a lattice for new tissue ingrowth.<sup>16</sup>
- As the scaffold bioabsorbs, the new ingrown tissue provides strength to the repair site.<sup>19</sup>
- By 52 weeks the new ingrown tissue is approximately 2.4 mm thick and provides a majority of strength to the repair site.<sup>2</sup>

By providing a lattice for tissue regeneration, GalaFLEX® encourages cells to migrate into its pores, allowing stronger, organized collagen to build and healthy blood vessels to form.<sup>1,16</sup>

G = GalaFLEX® scaffold • Human Tissue Specimen • Images shown at 100x magnification



Arrows denote new collagen formation

Arrows denote new blood vessels

**By 6 Weeks:**  
New tissue with abundant mature collagen (as indicated by positive type I collagen staining) and vascularization (as shown by positive CD31 and smooth muscle actin stains) has quickly integrated into the scaffold.<sup>1</sup>

**By 7 Months:**  
A fully integrated tissue plane of primarily type I collagen throughout the scaffold indicates collagen maturation and soft tissue regeneration (minimal inflammatory response with no evidence of encapsulation).<sup>1</sup>

## History of P4HB Products

1980s

Researchers at MIT developed a recombinant system to produce Polyhydroxyalkanoates (PHAs) in microorganisms.

1990s

Researchers at Metabolix further developed recombinant systems for the industrial production of PHAs. In 1998, Tepha, Inc. was incorporated to pursue the medical applications of PHAs.

2007 / 2008

The first P4HB medical devices: TephaFLEX® Suture & Mesh received FDA clearance.

2009 / 2010

Tepha partnered with B. Braun Medical who received the CE Mark for the P4HB device: MonoMax® Suture. MonoMax Suture was the first commercial launch of a P4HB device in Europe and the US.

2011

TephaFLEX Mesh received FDA clearance for soft tissue reinforcement in Plastic Surgery and was first used for Plastic Surgery. Tepha partnered with Tornier® and commercially launched: BioFiber™ for soft tissue reinforcement in the US.

2012 / 2013

Tepha partnered with Bard/Davol® to commercially launch the P4HB device: Phasix™ mesh for Hernia Repair in the US. Galatea Surgical, Inc.® became a wholly owned subsidiary of Tepha, Inc. to focus on plastic and reconstructive surgery.

2014 / 2015

Tepha P4HB devices achieved milestone of treating 1 million patients globally, with over 1,000 aesthetic plastic surgery patients. Galatea Surgical received CE Mark for use of GalaFLEX scaffold in breast surgery.

2016 / 2017

Galatea Surgical received FDA Clearance as the first and only 3-Dimensional scaffolds designed for plastic and reconstructive surgery.