## GalaFLEX® Scaffold

### Instructions for Use

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**ENGLISH****DESCRIPTION**

GalaFLEX® scaffold is a single layer, bioresorbable surgical mesh made from poly-4-hydroxybutyrate (P4HB). P4HB is produced from a naturally occurring monomer and is processed into monofilament fibres and knitted into a surgical scaffold. The structure of the scaffold consists of a knit pattern that mitigates the potential for unravelling when the scaffold is cut. The scaffold was fabricated specifically to provide strength throughout the critical wound healing period in soft tissue repair and reinforcement indications.

The design of the scaffold allows for an open pore pattern that encourages the proliferation of healthy tissue at the repair site, after which the product completely bioresorbs. Although the scaffold loses strength with time, its porous construction was designed to allow native tissue ingrowth and gradual transfer of load from the scaffold to the tissue.

**INTENDED USE**

GalaFLEX® scaffold is intended for use, as an adjunct to sutures, for the reinforcement and repair of soft tissue where weakness exists and where the addition of a reinforcing material is needed to obtain the desired surgical result in patients undergoing breast surgery.

The GalaFLEX® scaffold is designed to be used in patients undergoing soft tissue repair and reinforcement in medically necessary breast surgery procedures where the existing soft tissue is deficient to support the surgical repair. Examples of such breast surgery applications include reduction mammoplasty and breast revision surgery to correct a medical condition. GalaFLEX® scaffold may also be used in cosmetic breast procedures.

**CONTRAINdications**

None known.

**WARNINGS**

1. Device manufacture involves exposure to tetracycline hydrochloride and kanamycin sulfate. The safety and product use for patients with hypersensitivities to these antibiotics is unknown.
2. The safety and effectiveness of GalaFLEX® scaffold in neural tissue and in cardiovascular tissue has not been established.
3. The safety and effectiveness of GalaFLEX® scaffold in paediatric use has not been established.
4. Placement of the scaffold in direct contact with bowel or viscera is not recommended.
5. If an infection develops, treat the infection aggressively. An unresolved infection may require removal of the scaffold.
6. Because GalaFLEX® scaffold is fully bioresorbable, it should not be used in repairs where permanent support from the scaffold is required.
7. GalaFLEX® scaffold is supplied sterile. Inspect the device and packaging prior to use to be sure they are intact and undamaged.

8. GalaFLEX® is for single use only. Do not re-sterilise or re-use any portion of GalaFLEX® scaffold. Never re-use a scaffold, even if it seems to be in perfect condition, in order to prevent risk of cross-contamination or risk of reduced performance. Do not re-sterilise the scaffold. Once the package has been opened, any unused scaffold should be discarded.
9. Unused scaffold must be discarded according to the institution's procedures for handling of biohazardous materials.
10. Any decision to remove the scaffold should take into account potential risks associated with a second surgical procedure. This may include difficulty to remove the scaffold due to the ingrowth of tissue or degradation of scaffold. Scaffold removal should be followed by adequate post-operative management.

## **PRECAUTIONS**

Only doctors qualified in the appropriate surgical techniques should use this device. Users should be familiar with surgical procedures and techniques, including strength requirements and scaffold size choices. Improper selection, placement, positioning and fixation of GalaFLEX® can cause subsequent undesirable results. Patient should be provided with instructions regarding post-operative care, i.e. lifting, hygiene, activity limitations and any other specific patient/procedure requirements. GalaFLEX® scaffold is bioresorbable and therefore any initial palpability of the scaffold decreases over time.

## **ACTIONS**

GalaFLEX® scaffold degrades through a process of hydrolysis and hydrolytic enzymatic digestion process.

It has been developed to minimise the variability of resorption rate and strength and provide support throughout the expected period of healing.

Pre-clinical implantation studies indicate that GalaFLEX® scaffold retains approximately 70% of its strength at 12 weeks. Absorption of the scaffold material will be essentially complete within 18-24 months.

## **ADVERSE REACTIONS**

In pre-clinical testing, GalaFLEX® scaffold elicited a minimal tissue reaction characteristic of foreign body response to a substance. The tissue reaction resolved as the scaffold was resorbed. Possible complications of using GalaFLEX® scaffold include infection, seroma, pain or swelling, scaffold migration, wound dehiscence, haemorrhage, adhesions, haematoma, inflammation, extrusion and recurrence of the soft tissue defect.

## **IMAGING COMPATIBILITY**

GalaFLEX® scaffold is made from fully resorbable P4HB and is temporary in nature. This material is non-conducting, non-metallic and non-magnetic. Therefore, in accordance with the definition stated in ASTM F-2503-13, Standard Practice for Marking Devices and Other Items for Safety in the Magnetic Resonance Environment, the device is determined to be "MR Safe - an item that poses no known hazards in all MR environments."

## **DIRECTIONS FOR USE**

1. Prior to implanting GalaFLEX® scaffold, complete repair using appropriate suture technique.
2. GalaFLEX® scaffold should be prepared for placement using standard surgical preparation technique.
3. Using aseptic technique, GalaFLEX® scaffold may be cut to the shape or size desired for a specific application.
4. To ensure the scaffold remains flat, minimise pulling and stretching of the scaffold prior to implantation.
5. Transfer the scaffold to the surgical site. Implant the scaffold so that it lays flat against the tissue and the edge of the scaffold extends beyond the margins of the defect. Suture the scaffold into place, avoiding excessive tension and without the expectation of stretch.
6. It is recommended that suture fixation be placed 6 mm to 12 mm apart at a distance approximately 6 mm from the edge of the scaffold using an interrupted suturing technique.
7. If using a resorbable suture to fixate the scaffold, please ensure the strength retention of the suture is sufficient to secure the scaffold during the healing period.
8. Close the incision site using standard surgical technique.
9. Discard any unused portions of the scaffold per your institution's procedures for biohazardous materials.

## **HOW SUPPLIED**

GalaFLEX® scaffold is available in single packets as a sterile, undyed scaffold in single sheet sizes of varying widths and lengths.

## **STERILISATION**

GalaFLEX® scaffold is sterilised using Ethylene Oxide. GalaFLEX® scaffold is supplied sterile and is for SINGLE USE ONLY.

DO NOT clean, re-sterilise or re-use as this may damage or compromise performance of devices and may expose patient to risk of transmitting infectious disease.

## **STORAGE**

Store at room temperature, 15° to 25°C. Avoid prolonged exposure to elevated temperatures.

## SYMBOLS KEY

|                   |  |
|-------------------|--|
| <b>REF</b>        | Catalogue Number                                   |
| <b>LOT</b>        | Batch Code   |
|                   | Use By – year, month and day                       |
|                   | Do Not Re-use                                      |
|                   | Do not use if package is damaged                   |
|                   | Temperature limitation                             |
| <b>Rx Only</b>    | Rx ONLY Prescription Use Only                      |
| <b>STERILE EO</b> | Sterilised by Ethylene Oxide                       |
|                   | Do Not Re-sterilise                                |
|                   | CAUTION - Consult Instructions for Use             |
|                   | Manufacturer                                       |
| <b>EC REP</b>     | European Union Authorised Representative           |
| <b>CE 2797</b>    | CE mark and identification number of Notified Body |
| <b>MR</b>         | MR Safe  |

## CONTACT INFORMATION

Manufacturer

Tepha, Inc.  
99 Hayden Avenue, Suite 360  
Lexington, MA 02421 USA  
Galatea Surgical, Inc. is a wholly owned  
subsidiary of Tepha, Inc.

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Galatea Surgical, Inc., tamamına Tepha Inc.'in